


## Product Datasheet

### Human CCK ELISA Kit (orb562418)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Human CCK ELISA Kit                    |
| <b>Reactivity</b>        | Human                                  |
| <b>Range</b>             | 15.6-1000pg/ml                         |
| <b>Target</b>            | CCK                                    |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Uniprot ID</b>        | <b>P06307</b>                          |
| <b>Sensitivity</b>       | 9.375pg/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Human  
CCK(Cholecystokinin)  
ELISA Kit