


## Product Datasheet

### Human Mucin 4 ELISA Kit (orb562546)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Human Mucin 4 ELISA Kit |
| <b>Reactivity</b>      | Human                   |
| <b>Range</b>           | 0.313-20ng/ml           |
| <b>Target</b>          | MUC4                    |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>Q99102</b>           |
| <b>Sensitivity</b>     | 0.188ng/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Human  
MUC4(Mucin-  
4) ELISA Kit