


## Product Datasheet

### Human Syntaxin 1A ELISA Kit (orb562678)

|                        |                             |
|------------------------|-----------------------------|
| <b>Description</b>     | Human Syntaxin 1A ELISA Kit |
| <b>Reactivity</b>      | Human                       |
| <b>Range</b>           | 0.156-10ng/ml               |
| <b>Target</b>          | STX1A                       |
| <b>Storage</b>         | 2-8 °C for 6 months         |
| <b>Note</b>            | For research use only       |
| <b>Uniprot ID</b>      | <b>Q16623</b>               |
| <b>Sensitivity</b>     | 0.094ng/ml                  |
| <b>Expiration Date</b> | Please enquire.             |



Human  
STX1A(Syntaxin-  
1A) ELISA Kit