


Product Datasheet

Human B3GALT1 ELISA Kit (orb562741)

Description	Human B3GALT1 ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	B3GALT1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9Y5Z6
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Human
B3GALT1(Beta-1,
3-
galactosyltransf...