


Product Datasheet

Human B3GALT2 ELISA Kit (orb562743)

Description	Human B3GALT2 ELISA Kit
Reactivity	Human
Range	0.312-20ng/ml
Target	B3GALT2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	O43825
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
B3GALT2(Beta-1,
3-
galactosyltransf...