


## Product Datasheet

### Human C19orf80 ELISA Kit (orb562745)

|                        |                          |
|------------------------|--------------------------|
| <b>Description</b>     | Human C19orf80 ELISA Kit |
| <b>Reactivity</b>      | Human                    |
| <b>Range</b>           | 46.875-3000pg/ml         |
| <b>Target</b>          | C19orf80                 |
| <b>Storage</b>         | 2-8 °C for 6 months      |
| <b>Note</b>            | For research use only    |
| <b>Uniprot ID</b>      | <b>Q6UXH0</b>            |
| <b>Sensitivity</b>     | 28.125pg/ml              |
| <b>Expiration Date</b> | Please enquire.          |



Human  
C19orf80(Betatrophin)  
ELISA Kit