


Product Datasheet

Human AP1M1 ELISA Kit (orb562839)

Description	Human AP1M1 ELISA Kit
Reactivity	Human
Range	0.312-20ng/ml
Target	AP1M1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9BX55
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
AP1M1(AP-
1 complex
subunit
mu-1) E...