


Product Datasheet

Human RRM2B ELISA Kit (orb562976)

Description	Human RRM2B (Ribonucleoside-diphosphate reductase subunit M2 B)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	RRM2B
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7LG56
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human RRM2B
(Ribonucleoside-
diphosphate ...