


## Product Datasheet

### Human LOXL2 ELISA Kit (orb563127)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Human LOXL2 ELISA Kit                  |
| <b>Reactivity</b>        | Human                                  |
| <b>Range</b>             | 37.5-2400pg/ml                         |
| <b>Target</b>            | LOXL2                                  |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Uniprot ID</b>        | <b>Q9Y4K0</b>                          |
| <b>Sensitivity</b>       | 22.5pg/ml                              |
| <b>Expiration Date</b>   | Please enquire.                        |



Human  
LOXL2(Lysyl  
oxidase  
homolog 2)  
ELI...