


Product Datasheet

Human PLEKHO1 ELISA Kit (orb563153)

Description	Human PLEKHO1 ELISA Kit
Reactivity	Human
Range	0.625-40ng/ml
Target	PLEKHO1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q53GL0
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Human
PLEKHO1(Pleckstrin
homology
domain...