


Product Datasheet

Human TTPAL ELISA Kit (orb563156)

Description	Human TTPAL ELISA Kit
Reactivity	Human
Range	1.56-100ng/ml
Target	TTPAL
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9BTX7
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
TTPAL(Alpha-
tocopherol
transfer pr...