


## Product Datasheet

### Human SLC25A4 ELISA Kit (orb563164)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Human SLC25A4 ELISA Kit |
| <b>Reactivity</b>      | Human                   |
| <b>Range</b>           | 0.156-10ng/ml           |
| <b>Target</b>          | SLC25A4                 |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>P12235</b>           |
| <b>Sensitivity</b>     | 0.094ng/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Human  
SLC25A4(ADP/ATP  
translocase 1)  
ELI...