


Product Datasheet

Human HLA-DRB5 ELISA Kit (orb563177)

Description	Human HLA-DRB5 ELISA Kit
Reactivity	Human
Range	1.56-100ng/ml
Target	HLA-DRB5
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q30154
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human HLA-DRB5(HLA class II histocompati...