

Product Datasheet

Zap-70 (phospho-Tyr315/Tyr319) antibody (orb5633)

Description	Zap-70 (phospho-Tyr315/Tyr319) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Zap-70 around the phosphorylation site of Tyr315/Tyr319 (SV(p-Y)ESP(p-Y)SD)
Target	Phospho-ZAP70 (Tyr315 + Tyr319)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68 kDa
Uniprot ID	P43403
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.