


## Product Datasheet

### Human SIK2 ELISA Kit (orb563657)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Human SIK2 ELISA Kit                   |
| <b>Reactivity</b>        | Human                                  |
| <b>Range</b>             | 0.312-20ng/ml                          |
| <b>Target</b>            | SIK2                                   |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Uniprot ID</b>        | <b>Q9H0K1</b>                          |
| <b>Sensitivity</b>       | 0.188ng/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Human  
SIK2(Serine/threonine-  
protein kina...