


Product Datasheet

Human TSPAN14 ELISA Kit (orb563731)

Description	Human TSPAN14 ELISA Kit
Reactivity	Human
Range	0.312-20ng/ml
Target	TSPAN14
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q8NG11
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
TSPAN14(Tetraspanin-
14) ELISA Kit