


Product Datasheet

Human MTPAP ELISA Kit (orb563742)

Description	Human MTPAP ELISA Kit
Reactivity	Human
Range	0.312-20ng/ml
Target	MTPAP
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9NVV4
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
MTPAP(Poly(A)
RNA
polymerase,
mito...