


Product Datasheet

Human C9orf140 ELISA Kit (orb563749)

Description	Human C9orf140 ELISA Kit
Reactivity	Human
Range	0.78-50ng/ml
Target	C9orf140
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q86UD0
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Human
C9orf140(Protein
C9orf140) ELISA
K...