


## Product Datasheet

### Human PLA2G2A ELISA Kit (orb563794)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Human PLA2G2A ELISA Kit |
| <b>Reactivity</b>      | Human                   |
| <b>Range</b>           | 0.938-60ng/ml           |
| <b>Target</b>          | PLA2G2A                 |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>P14555</b>           |
| <b>Sensitivity</b>     | 0.563ng/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Human  
PLA2G2A(Phospholipase  
A2, membrane...