


## Product Datasheet

### Human PDX1 ELISA Kit (orb563806)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human PDX1 ELISA Kit  |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 15.6-1000pg/ml        |
| <b>Target</b>          | PDX1                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P52945</b>         |
| <b>Sensitivity</b>     | 9.375pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
PDX1(Pancreas/duodenum  
homeobox pr...