


## Product Datasheet

### Human C4BPbeta ELISA Kit (orb563998)

|                        |                          |
|------------------------|--------------------------|
| <b>Description</b>     | Human C4BPbeta ELISA Kit |
| <b>Reactivity</b>      | Human                    |
| <b>Range</b>           | 0.469-30ng/ml            |
| <b>Target</b>          | C4BPβ                    |
| <b>Storage</b>         | 2-8 °C for 6 months      |
| <b>Note</b>            | For research use only    |
| <b>Uniprot ID</b>      | <b>P20851</b>            |
| <b>Sensitivity</b>     | 0.281ng/ml               |
| <b>Expiration Date</b> | Please enquire.          |



Human  
C4BPβ(C4  
Binding  
Protein  
Beta) EL...