


## Product Datasheet

### Human DBP ELISA Kit (orb564098)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Human DBP ELISA Kit     |
| <b>Reactivity</b>      | Human                   |
| <b>Range</b>           | 0.156-10ng/ml           |
| <b>Target</b>          | DBP/DNA Binding Protein |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>Q10586</b>           |
| <b>Sensitivity</b>     | 0.094ng/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Human  
DBP/DNA  
Binding  
Protein(DNA  
Bindin...