


Product Datasheet

Human FcgammaR3A ELISA Kit (orb564156)

Description	Human FcgammaR3A(Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
Reactivity	Human
Range	7.813-500ng/ml
Target	FcγR3A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P08637
Sensitivity	4.688ng/ml
Expiration Date	Please enquire.



Human
FcγR3A(Fc
Fragment of
IgG Low
Aff...