

Product Datasheet

Human ITG $\alpha V\beta 3$ (Integrin $\alpha V\beta 3$) ELISA Kit (orb564178)

Description Human ITG $\alpha V\beta 3$ (Integrin $\alpha V\beta 3$) ELISA Kit

Reactivity Human

Range 31.2-2000pg/ml

Target ITG $\alpha V\beta 3$


Storage 2-8 °C for 6 months

Note For research use only

Uniprot ID **P06756, P05106**

Sensitivity 18.75pg/ml

Expiration Date Please enquire.



Human ITG $\alpha V\beta 3$ (Integrin $\alpha V\beta 3$) EL...