

Product Datasheet

AMY2A antibody (orb5642)

Description	AMY2A antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, WB
Immunogen	KLH conjugated synthetic peptide derived from human Pancreatic Amylase (421-511/511 aa)
Target	AMY2A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P04746
Expiration Date	12 months from date of receipt.