


Product Datasheet

Human GSK3alpha ELISA Kit (orb564234)

Description	Human GSK3alpha(Glycogen Synthase Kinase 3 Alpha) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	GSK3 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P49840
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
GSK3 α (Glycogen
Synthase Kinase
3 ...