



Product Datasheet

14-3-3 zeta antibody (orb5643)

Description	14-3-3 zeta antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human 14-3-3 protein zeta/delta around the phosphorylation site of Ser58 (RS(p-S)WR)
Target	phospho-14-3-3 protein zeta (Ser58)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27 kDa
Uniprot ID	P63104
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A549 Cell
US...



Immunohistochemical staining of human
br...