


Product Datasheet

Human MIP4alpha ELISA Kit (orb564324)

Description	Human MIP4alpha(Macrophage Inflammatory Protein 4 Alpha) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	MIP4 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y258
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
MIP4 α (Macrophage
Inflammatory
Pro...