


## Product Datasheet

### Human MUC3 ELISA Kit (orb564339)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human MUC3 ELISA Kit  |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 31.25-2000pg/ml       |
| <b>Target</b>          | MUC3                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q02505</b>         |
| <b>Sensitivity</b>     | 18.75pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
MUC3(Mucin  
3) ELISA Kit