


Product Datasheet

Human PGF2alpha ELISA Kit (orb564420)

Description	Human PGF2alpha(Prostaglandin F2alpha) ELISA Kit
Reactivity	Human
Range	7.813-500pg/ml
Target	PGF2 α
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Human
PGF2 α (Prostaglandin
F2 α) ELISA K...