


Product Datasheet

Human PKC ζ (Protein Kinase C Zeta) ELISA Kit (orb564430)

Description	Human PKC ζ (Protein Kinase C Zeta) ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	PKC ζ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q05513
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
PKC ζ (Protein
Kinase C
Zeta)
ELISA...