


Product Datasheet

Human PLCgamma2 ELISA Kit (orb564443)

Description	Human PLCgamma2(Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	PLCγ2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P16885
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
PLCγ2(Phospholipase
C Gamma 2, Ph...