


Product Datasheet

Human PPAR-alpha ELISA Kit (orb564455)

Description	Human PPAR-alpha(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PPAR- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q07869
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human PPAR- α (Peroxisome Proliferator Ac...