


Product Datasheet

Human PPAR-gamma ELISA Kit (orb564456)

Description	Human PPAR-gamma(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PPAR-γ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P37231
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PPAR-
γ(Peroxisome
Proliferator
Ac...