


Product Datasheet

Human PTF1alpha ELISA Kit (orb564479)

Description	Human PTF1alpha(Pancreas Specific Transcription Factor 1alpha) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PTF1 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7RTS3
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PTF1 α (Pancreas
Specific
Transcrip...