


## Product Datasheet

### Human RPS6Kbeta1 ELISA Kit (orb564501)

|                        |                            |
|------------------------|----------------------------|
| <b>Description</b>     | Human RPS6Kbeta1 ELISA Kit |
| <b>Reactivity</b>      | Human                      |
| <b>Range</b>           | 0.156-10ng/ml              |
| <b>Target</b>          | RPS6Kβ1                    |
| <b>Storage</b>         | 2-8 °C for 6 months        |
| <b>Note</b>            | For research use only      |
| <b>Uniprot ID</b>      | <b>P23443</b>              |
| <b>Sensitivity</b>     | 0.094ng/ml                 |
| <b>Expiration Date</b> | Please enquire.            |



Human  
RPS6Kβ1(Ribosomal  
Protein S6 Kina...