


Product Datasheet

Human TMSbeta4 ELISA Kit (orb564569)

Description	Human TMSbeta4 ELISA Kit
Reactivity	Human
Range	31.25-2000ng/ml
Target	TMSβ4
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P62328
Sensitivity	18.75ng/ml
Expiration Date	Please enquire.



Human
TMSβ4(Thymosin
Beta 4) ELISA Kit