


## Product Datasheet

### Human NXPE4 ELISA Kit (orb564778)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human NXPE4 ELISA Kit |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 0.156-10ng/ml         |
| <b>Target</b>          | NXPE4                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q6UWF7</b>         |
| <b>Sensitivity</b>     | 0.094ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
NXPE4(NXPE  
family  
member 4)  
ELISA ...