

Product Datasheet

Human anti-HSPA5 ELISA Kit (orb564926)

Description Human anti-HSPA5(anti-Heat shock protein family A member 5)ELISA KIT

Reactivity Human

Range 0.391-25ng/ml

Target anti-HSPA5

Storage 2-8 °C for 6 months


Note For research use only

Application notes indirect elisa

Uniprot ID **P11021**

Sensitivity 0.234ng/ml

Expiration Date Please enquire.



Human anti-HSPA5(anti-Heat shock protein...