


Product Datasheet

Human CHCHD2 ELISA Kit (orb564959)

Description	Human CHCHD2 (Coiled-coil-helix-coiled-coil-helix domain-containing protein 2) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CHCHD2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y6H1
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CHCHD2(Coiled-
coil-helix-coiled-
co...