


## Product Datasheet

### Human ERBB4 ELISA Kit (orb564974)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human ERBB4 ELISA Kit |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 0.313-20ng/ml         |
| <b>Target</b>          | ERBB4                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q15303</b>         |
| <b>Sensitivity</b>     | 0.188ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
ERBB4  
(Receptor  
tyrosine-  
protein  
k...