


Product Datasheet

Human IGF2BP1 ELISA Kit (orb565011)

Description	Human IGF2BP1(Insulin-like growth factor 2 mRNA-binding protein 1)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	IGF2BP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9NZI8
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
IGF2BP1(Insulin-
like growth
factor...