


Product Datasheet

Hamster TNF-alpha ELISA Kit (orb565036)

Description	Hamster TNF-alpha(Tumor Necrosis Factor-alpha) ELISA Kit
Reactivity	Hamster
Range	3.12-200pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A0A1U7RE37
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.



Hamster
TNF-
 α (Tumor
Necrosis
Factor- α)...