


## Product Datasheet

### Hamster IL-10 ELISA Kit (orb565040)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Hamster IL-10 ELISA Kit |
| <b>Reactivity</b>      | Hamster                 |
| <b>Range</b>           | 6.25-400pg/ml           |
| <b>Target</b>          | IL-10                   |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>G3HB01</b>           |
| <b>Sensitivity</b>     | 3.75pg/ml               |
| <b>Expiration Date</b> | Please enquire.         |



Hamster IL-10(Interleukin-10) ELISA Kit