


## Product Datasheet

### Mouse TIMP-4 ELISA Kit (orb565168)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Mouse TIMP-4 ELISA Kit |
| <b>Reactivity</b>      | Mouse                  |
| <b>Range</b>           | 62.5-4000pg/ml         |
| <b>Target</b>          | TIMP-4                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>Q9JHB3</b>          |
| <b>Sensitivity</b>     | 37.5pg/ml              |
| <b>Expiration Date</b> | Please enquire.        |



Mouse  
TIMP-  
4(Tissue  
Inhibitors  
of  
Metall...