


Product Datasheet

Mouse TNFalpha ELISA Kit (orb565170)

Description	Mouse TNFalpha(Tumor Necrosis Factor Alpha) ELISA Kit
Reactivity	Mouse
Range	3.906-250pg/ml
Target	TNF α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P06804
Sensitivity	2.344pg/ml
Expiration Date	Please enquire.



Mouse
TNF α (Tumor
Necrosis
Factor
Alpha)...