


Product Datasheet

Mouse Neuropeptide Y ELISA Kit (orb565216)

Description	Mouse Npy(neuropeptide Y) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Npy
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Sandwich ELISA, Double Antibody
Uniprot ID	P57774
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
Npy(Pro-
neuropeptide
Y) ELISA Kit