


## Product Datasheet

### Mouse Aph1a ELISA Kit (orb565295)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse Aph1a ELISA Kit |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 1.56-100ng/ml         |
| <b>Target</b>          | Aph1a                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q8BVF7</b>         |
| <b>Sensitivity</b>     | 0.938ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
Aph1a(Gamma-  
secretase  
subunit APH-...