


Product Datasheet

Mouse Gabra2 ELISA Kit (orb565302)

Description	Mouse Gabra2 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Gabra2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P26048
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Gabra2(Gamma-aminobutyric acid rec...