


Product Datasheet

Mouse SOD1 ELISA Kit (orb565320)

Description	Mouse SOD1 ELISA Kit
Reactivity	Mouse
Range	0.78-50ng/ml
Target	SOD1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P08228
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Mouse
SOD1(Superoxide
Dismutase [Cu-
Zn])...